



**VEGA-6 /
150 mA
HIGH FREQUENCY
FIXED X-RAY MACHINE**

About Innovation Meditech (P) Ltd.

Innovation Meditech (P) Ltd. Is established by team of professional with decades of experience in healthcare technology, with a vision to offer products and services of global standards. The manufacturing facility is spread over 25000 sq. Feet with stringent quality control systems. All Products Manufactured are in compliance with the standards provided by Bureau of Indian Standards (BIS) for Mechanical and Electrical Safety and Atomic Energy Regulatory Board (AERB). Innovation Meditech (P) Ltd. Has nationwide presence with network of sales offices, resident sales representatives and dealers.

Innovation Service Advantage

Innovation Meditech (P) Ltd. is AERB approved Services agency having nationwide presence of qualified and trained service engineers to take care of installation, maintenance, repair and quality assurance test of all type of Radiology equipments.

OTHER PRODUCTS

- ➔ 100 mA Mobile X-Ray
- ➔ 100 mA Fix X-ray
- ➔ 300 mA X-ray
- ➔ 500 mA X-ray
- ➔ OPG Machine
- ➔ 4.0kW HF Mobile X-ray Machine
- ➔ 7.0kW HF Mobile X-Ray Machine
- ➔ 15.0 HF X-Ray Machine
- ➔ 32kW HF X-Ray Machine
- ➔ 40kW HF X-Ray Machines

SERVICES AVAILABLE

- ➔ AERB and NABL accreditation Quality Assurance Test for all Type of X-Ray & CT Scan Machine.



INNOVATION MEDITECH PVT. LTD.

Toll Free No.: 1800-2026-276

Phone No.: 0120-499-6085

Address : B-6, Sector-88, Noida, Gautam Budh Nagar - 201305, UP, (INDIA)

E-Mail : info@innovationmeditech.com

Website : www.innovationmeditech.com



BRANCH OFFICE & RESIDENT

JAIPUR, LUCKNOW, PATNA, HYDERABAD, HISAR, KOLKATA



www.innovationmeditech.com

Certified and Approved By:





TECHNICAL SPECIFICATION

S. No.	PARAMETER	SPECIFIED VALUE
1	Mains Voltage	230VAC ±10%
2	Frequency	50 Hz
3	Power Output	6 kW
4	Generator Frequency	40 kHz
5	mA Range	150mA
6	KV Range	120kV
7	kV Accuracy	±5kV
8	mAs Accuracy	±10%
9	X-ray Tube	Stationary Anode
10	Target Angle	15°
11	Focal Spot	2.0 mm
12	Total Filtration	2.5 mmAl Eqv.
13	Inherent Tube Filtration	Min. 0.65 mm Al / 75 kV
14	Collimator	Manual
15	Exposure Release Switch	2 Steps
16	Parameter Selection	kV & mAs
17	Mode	APR (Anatomical Programmed Radiography) & Manual
18	Conforms to Standadrds	IS:7620 Part-1
		AERB Type Approved
19	AERB Approval	Approved

FEATURES

- POWERED WITH HUMAN MACHINE INTERFACE
- LIGHT WEIGHT
- HF SMART GENERATOR
- SLEEK DESIGN
- AERB APPROVED



- * SAFETY STANDARD BIS : IS7620 : PART-1
- * REGULATORY REQUIREMENT : AERB TYPE APPROVAL

IMPROVE IMAGE QUALITY WITH A HIGH-PERFORMANCE CONTROL PANEL

